BIOMED - TP 1

Detection and tracking of spots

Marie Diez

1 Detect and track spots

Protocole to detect and track spots :

- Clean the image
 We can use a threshold plugin to keep the white spots like the *threshold* plugin using values 450 600.
- Use spot tracking
 First we have to check the swimming pool option for the tracking next, i also used scale 2 to detect
 object size around 3 pixels and i check the bright spot over black background option. Then after
 that the detection is done we have to track the spot, i use an expected track length of 10 because
 of my computer performances.

I use a cropped version of the image $Tracking_FITC.tif$ because of my computer performances, without that the detection crash every time.

This is some results :




Figure 1: Spots detection after thresholding




Figure 2: Spots tracking




Figure 3: Track lenght




Figure 4: Average length computation

We can also have lot more informations with all the track processor plugins :




Figure 5: Tracking performance computation




Figure 6: ROI Statistics




Figure 7: Motion profiler